

Or, as the person you really are?



Oxy-View is a product of Oxy-View Inc., a respected maker of proven medical devices which help people needing continuous oxygen therapy live better lives.

OXY-VIEW INC.

109 Inverness Drive East
Englewood, Colorado 80112-5105

Patients, eye care professionals, physicians, respiratory therapists and other interested parties are invited to contact us toll free at 877-699-8439 or online at cservice@oxyview.com and www.oxyview.com.

OXY-VIEW EYEWEAR ORDERING INFORMATION:

The first six digits of the order number indicate the frame style (Half Rim "HR" or Full Rim "FR") and the lens size (51/17 or 53/19). The Last four digits indicate the color of the frame. For example, HR5117/3116 is a Half-Rim frame, 51/17 lens size, color rose.

Order No.	Description	"B"
	<i>Oval, Half-Rim</i>	
HR5117/3116	Rose	36 mm
HR5117/3062	Metalic Light Blue	36 mm
HR5117/3081	Metalic Red	36 mm
HR5117/3059	Shaded Dark Blue	36 mm
HR5117/3008	Shaded Gold	36 mm
	<i>Retangular, Half-Rim</i>	
HR5319/3008	Shaded Gold	38 mm
HR5319/3068	Metalic Silver	38 mm
HR5319/3011	Shaded Dark Gunmetal	38 mm
HR5319/3091	Metalic Brown	38 mm
	<i>Retangular, Full-Rim</i>	
FR5319/3068	Metalic Silver	38 mm
FR5319/3091	Metalic Brown	38 mm

U.S. Patent 6,886,562

Do you want to
be seen as a person
on oxygen therapy?



Oxy-View eyewear lets people see the real you.

You need oxygen therapy. But you also need for people to still see you, and not just oxygen tubes. Now there's an answer—Oxy-View eyewear. These attractive eyeglass frames not only hold your prescription lenses, they also deliver the oxygen you need in a way that others will hardly notice. They'll see only your smiling face.

Oxygen therapy without the nasal cannula.


Let's face it – it's no fun to wear a nasal cannula with those tubes across your face and over your ears. But now you can wear Oxy-View eyeglasses instead. They can effectively deliver oxygen flows of up to five liters per minute, as prescribed by your physician. And they work with all types of oxygen delivery systems, including liquid oxygen, gas cylinders, concentrators, and most oxygen conserving devices.

Better comfort and a better look.


Oxy-View eyeglasses help you see and feel better. The frames come in a variety of styles and colors, and are perfect for most prescription and sunglass lenses. They offer you an attractive way of reducing the discomfort that often occurs from wearing a nasal cannula, as well as a boost to your spirits by disguising those unsightly tubes.

Go out without standing out.

When you're wearing Oxy-View eyeglasses, you'll feel better about getting out and doing the things you love. People will no longer see someone on oxygen, but rather the vibrant person you really are. And you'll see them clearly, too, with your Oxy-View eyeglasses.



A simple idea that works beautifully.



Choose from several distinctive frames.

Oxy-View eyeglass frames are available in several fashionable colors. Please call toll free 1-877-699-8439 for available colors. Whichever color you choose, you'll find them an attractive and comfortable way to get the oxygen your doctor has prescribed. Always consult your doctor before changing your oxygen delivery system.



Oval, Half-Rim, 51/17



Rectangular, Half-Rim, 53/19

Available Colors



Actual colors may vary.

See Rear Panel For Ordering Information